

SurgyCut™

NMSDC Certified • FDA Registered

A Reposable Laparoscopic Scissor

The SurgyCut reposable scissor is configured as a reusable handle and disposable shaft. Each of the three available SurgyCut scissor shafts connect seamlessly to our durable handle using a ball and socket connection, delivering a peak performance each time.

Available in both 35 cm and 45 cm shaft lengths.



Regular Metz Scissor Shaft

790-200-535: Reposable Scissor Shaft 5mm/35cm
790-200-545: Reposable Scissor Shaft 5mm/45cm



Mini Metz Scissor Shaft

790-202-535: Reposable Mini Scissor Shaft 5mm/35cm
790-202-545: Reposable Mini Scissor Shaft 5mm/45cm



Hook Scissor Shaft

790-204-535: Reposable Hook Scissor Shaft 5mm/35cm
790-204-545: Reposable Hook Scissor Shaft 5mm/45cm



Rotation Knob

Easy and smooth rotation with the index finger

Ergonomic Handle Design

Providing a comfortable grip when handling and cutting



U.S. SURGITECH
Your Source For Medical Products