


SurgyNeedle



The only disposable & clear 5mm aspiration needle!

Aspiration Solution: The SurgyNeedle

The SurgyNeedle is a patent pending aspirating needle used during laparoscopic and robotic procedures. Pulled for each case, but only used when needed, the SurgyNeedle comes sterile packed eliminating dullness, cleaning and sterilization. Our single use product minimizes the risk of contamination in the OR.



Sharp, Beveled Tip
Enables controlled cyst & organ penetration

Clear & Transparent Shaft
Provides greater fluid visualization

Semi-Rigid Shaft
Allows for more maneuverability

Luer-Lock Connection
Provides safe & simple insertion through a trocar

PRODUCT CODES

- 780-411: 15-gauge, 40 cm length, 5 mm
- 780-412: 21-gauge, 40 cm length, 5 mm
- 780-413: 18-gauge, 40 cm length, 5 mm
- 780-414: 23-gauge, 40 cm length, 5 mm
- 780-415: 23-gauge, 15 cm length, 5 mm