




Retrieval has never
been so easy.



Specimen Retrieval Bag

Laparoscopic Solutions: Specimen Retrieval Bag

The GeniStrong™ Bag is an easy-to-use retrieval system designed to temporarily contain specimen and facilitate its removal from the patient during laparoscopic surgery. Designed to minimize contamination of the abdominal cavity, the Specimen Retrieval Bag includes the following distinctive features:



Specimen Retrieval Bag

Description	Size (mouth x length)	Qty/Box	Volume	Part Number
Medium Specimen Retrieval Bag (10mm shaft)	7.3 cm x 17.0 cm	10	240 mL	550-000-003
Large Specimen Retrieval Bag (12mm shaft)	11 cm x 22 cm	10	750mL	550-000-005
Extra Large Specimen Retrieval Bag (15mm shaft)	14 cm x 27 cm	10	2000 mL	550-000-008