


Performance you expect,
Cost savings you need.

LiNA PowerBlade Plus™

Advanced Bipolar Sealer/Dissector



LiNA PowerBlade Plus™

Confidence, during the procedure and after.

PowerBlade Plus™ is capable of sealing vessels up to 7mm with seal strength in excess of 3x systolic pressure and minimal thermal spread to adjacent tissue.¹

Tip-first jaw closure facilitates grasping of even the smallest tissue bundles.

Low profile jaw with textured surface for delicate dissection of tissue planes.



Control, without the complexity.

We've intuitively designed the PowerBlade Plus™ with Clamping, Sealing, and Cutting as independent functions for improved safety and control in the hands of every surgeon.



Compatibility, with the generators you already own.


PowerBlade Plus™ is designed to work with nearly all standard bipolar ESU's so you won't need to clutter up the OR with another expensive generator.



How does PowerBlade Plus™ measure up to the competition?



Jaw and Cut Length Comparison¹
Actual Size



Ordering Information:

Product code	Description	Unit of measure
5S-SJ-370-YH	LiNA PowerBlade Plus™ - 5mm endoscopic bipolar coagulation and cutting forceps - Sterile, Single use.	Box of 6 ea

*Please refer to the IFU for Intended Use and Contraindications.

Distributor Contact:

LiNA Medical USA

1856 Corporate Drive, Suite 135 | Norcross, GA 30093 | Tel: (855) 546-2633 | Fax: (678) 380-0500 | orders@linamed.com | www.linamed.com
Manufactured by LiNA Medical Aps, Formervangen 5, 2600 Glostrup, Denmark

©2015 LiNA Medical USA Inc. LiNA PowerBlade Plus 7/2015 USA.. 1. Bench study, data on file.

LigaSure and Sonicision are Trademarks of Covidien. Enseal and Harmonic ACE are Registered Trademarks of Ethicon-Endo Surgery. PKS and Thunderbeat are Trademarks of Olympus.

