




# ElektroBlue™ INSTRUMENTS





# ElektroBlue™

Over a decade of design input and services has led us to develop a complete line of non conductive instrumentation called Elektroblue.

These instruments are manufactured in the United States with imported medical grade German forged stainless steel blanks. This gives us the capability of customizing any instrumentation with insulation and fast turnaround times. Striving for constant improvement, U.S. Surgitech only uses FDA approved nylon 11 powder, while implementing processes to become ISO 13485 compliant, and following strict guidelines including but not limited to batch hipot testing.

Each and every instrument is laser marked with a lot# to allow for traceability and a one year warranty free of manufacture defects.



# FORCEPS



Elektroblue Leep forceps range in a wide array of sizes and tips depending on your preference and use. The various tip choices with customizable lengths aid in grasping and removing tissue securely.

## **ADSON DRESSING FORCEPS**

750-101 Serrated 12.1cm/4.75in

## **ALLIS TISSUE FORCEPS**

750-102 5X6 19cm/7.5in

750-103 5X6 24.1cm/9.5in

## **BABCOCK FORCEPS**

750-104 20cm/8in

750-105 24.1cm/9.5in

## **BOZEMAN FORCEPS**

750-106 Straight 26cm/10.25in

750-107 S-Curved 26cm/10.25in

750-108 S-Curved One LG Ring 26cm/10.25in

## **CAMPION FORCEPS**

750-109 23cm/9in

## **CLEMETSON FORCEPS**

750-110 NO Ratchet 24.1cm/9.5in

## **DEBAKEY TISSUE FORCEPS**

750-111 1.5 mm jaws: 15.2cm/6in

750-112 1.5 mm jaws: 19.7cm/7.75in

750-113 1.5 mm jaws: 24.1cm/9.5in

750-114 1.5 mm jaws: 30.5cm/12in

750-115 2.0mm Jaws: 15.2cm/6in

750-116 2.0mm Jaws: 20cm/7.75in

750-117 2.0mm Jaws: 24cm/9.5in

750-118 2.0mm Jaws: 30.5cm/12in

750-119 2.0mm Jaws: 35.6cm/14in

750-120 2.7mm Jaws: 15cm/6in

750-121 2.7mm Jaws: 20cm/7.75in

750-122 2.7mm Jaws: 24cm/9.5in

750-123 2.7mm Jaws: 30.5cm/12in

750-124 2.7mm Jaws: 35cm/14in

750-125 2.7mm Jaws: 41cm/16in

750-126 2.7mm Jaws: 45.5cm/18in

## **DRESSING FORCEPS**

750-127 Serrated Exposed Tips 20cm/8in

750-128 Serrated Coated Tips 20cm/8in

750-129 Serrated Exposed Tips 25.5cm/10in

750-130 Serrated Coated Tips 25.5cm/10in

750-131 Serrated Exposed Serrated Tip 30.5cm/12in

750-132 Serrated Coated Serrated Tip 30.5cm/12in

750-133 Serrated Coated Tips 35.5cm/14in

750-134 Serrated Exposed Tips 35.5cm/14in

## **IRIS TISSUE FORCEPS**

750-135 1/2 Curved 10cm/4in

750-136 Straight 10cm/4in

## **LAUFE UTERINE FORCEPS**

750-137 Laufe Uterine Forceps 20.3cm/8in

## **POTTS SMITH FORCEPS**

750-138 1x2 Teeth 2.0mm Jaw 18cm/7in

750-139 1x2 Teeth 2.0mm Jaw 21cm/8.25in

750-140 1x2 Teeth 2.0mm Jaw 24cm/9.5in

750-141 2.0mm Jaws 18cm/7in

750-142 2.0mm Jaw 21cm/8.25in

750-143 2.0mm Jaw 24cm/9.5in

## **RING FORCEPS**

750-144 24cm/9.7in

## **SINGLEY FORCEPS**

750-145 23cm/9in

## **TISSUE FORCEPS**

750-146 1x2 Teeth Exposed Tips 17cm/6.75in

750-147 1x2 Teeth Coated Tips 17cm/6.75in

750-148 1x2 Teeth Exposed Tips 20cm/8in

750-149 1x2 Teeth Coated Tips 20cm/8in

750-150 1x2 Teeth Exposed Tips 25.5cm/10in

750-151 1x2 Teeth Coated Tips 25.5cm/10in

## **TUTTLE FORCEPS**

750-152 18cm/7in

# HOOKS



Elektroblue Hooks have strong sharp tips w/various angles for optimum grip

## **EMMET HOOK**

750-201 Acute Angled 23cm/9in

750-202 Angled 23cm/9in

750-203 Double Angled 23cm/9in

750-204 Fully Curved SM 23cm/9in

750-205 Fully Curved LG 23cm/9in

## **IRIS HOOK**

750-206 23cm/9in

## **TWO PRONG HOOK**

750-208 Round handle 23cm/9in

## **THREE PRONG HOOK**

750-207 Round handle 23cm/9in

# RULERS



Elektroblue Leep Rulers are designed to offer the most precise measurements w/easy to read mm markings for depth

## **DEPTH RULER**

750-301 240mm/9.5in

## **WIDTH RULER**

750-302 210mm/8.25in



# TENACULUM



Elektroblue tenaculum is a prominent instrument that offers optimal stabilization of the surgical site.

## **SCHROEDER TENACULUM**

750-401 Exposed Tips 25cm/10in

750-402 Coated Tips 25cm/10in

## **BRAUN TENACULUM LITE TOUCH**

750-403 Exposed Tips 25cm/10in

# ENDOCERVICAL SPECULA



Elektroblue Endocervical Specula are designed for physicians to have an accurate view of the endocervical canal

## **KOGAN ENDOCERVICAL SPECULUM**

750-501 W/Spring Wide 28cm/11in

750-502 W/Spring Narrow 28cm/11in

750-503 Standard Tips 28cm/11in

750-504 Narrow Tips 28cm/11in

750-505 Screw Style Standard Tips 28cm/11in

750-506 Screw Style Narrow Tips 28cm/11in



# SPECULA

A wide array of Elektrobue Leep Speculums in many different shapes and sizes allow for easy access and viewing of the surgical site, available with your choice of smoke evacuation options.

## AUVARD SPECULUM

- 750-601 LEEP / LLETZ 1.1 Kg/2.5lb 10cm x 4.5cm/ 4in x 1.75in
- 750-602 LEEP / LLETZ 1.1 Kg/ 2.5 lb 8.5cm x 4cm/3.5in x 1.5in

## CLINTON

- 750-605 LEEP / LLETZ Graves Speculum Extra Long W/DSE Tube 15cm x 4.4cm/6in x 1.5in
- 750-606 LEEP / LLETZ Pederson Speculum Extra Long 15cm x 2.5cm/6in x 1in
- 750-607 LEEP / LLETZ Pederson Speculum Extra Long W/Smoke Tube 15cm x 2.5cm/6in x 1in
- 750-608 LEEP / LLETZ Pederson Speculum Extra Long W/DSE Tube 15cm x 2.5cm/6in x 1in

## COLLIN SPECULUM

- 750-609 LEEP / LLETZ Sm 9cm x 2.5cm/3.5in x 1in
- 750-610 LEEP / LLETZ Sm W/Smoke Tube 9cm x 2.5cm/3.5in x 1in
- 750-611 LEEP / LLETZ Sm W/DSE Tube 9cm x 2.5cm/3.5in x 1in
- 750-612 LEEP / LLETZ Med 10cm x 4cm / 4in x 1.5in
- 750-613 LEEP / LLETZ Med W/Smoke Tube 10cm x 4cm /4in x 1.5in
- 750-614 LEEP / LLETZ Med W/DSE Tube 10cm x 4cm / 4in x 1.5in
- 750-615 LEEP / LLETZ Lrg 11cm x 3.5cm/ 4.25in x 1.5in
- 750-616 LEEP / LLETZ Lrg W/Smoke Tube 11cm x 3.5cm/ 4.25in x 1.5in
- 750-617 LEEP / LLETZ Lrg W/DSE Tube 11cm x 3.5cm/ 4.25in x 1.5in

## CUSCO

- 750-618 LEEP / LLETZ Sm 2.2cm x 9.2cm/ .875in x 3.625in
- 750-619 LEEP / LLETZ Sm W/Smoke Tube 2.2cm x 9.2cm/ .875in x 3.625in
- 750-620 LEEP / LLETZ Sm W/DSE Tube 2.2cm x 9.2cm/.875in x 3.625in
- 750-621 LEEP / LLETZ Med 2.5cm x 10.2cm/1inx 4in
- 750-622 LEEP / LLETZ Med W/Smoke Tube 2.5cm x 10.2cm/1in x 4in
- 750-623 LEEP / LLETZ Med W/DSE Tube 2.5cm x 10.2cm/1inx 4in
- 750-624 LEEP / LLETZ Lrg 2.7cm x 11.1cm/ 1.0625" x 4.375in
- 750-625 LEEP / LLETZ Lrg W/Smoke Tube 2.7cm x 11.1cm/ 1.0625in x 4.375in
- 750-626 LEEP / LLETZ Lrg W/DSE Tube 2.7cm x 11.1cm/ 1.0625in x 4.375in

## FREEWAY SPECULUM

- 750-627 LEEP / LLETZ Graves Speculum Open 90Deg
- 750-628 LEEP / LLETZ Speculum W/Smoke Tube 45Deg

## GRAVES SPECULUM

- 750-629 LEEP / LLETZ X Sm 7.5cmx1.5cm/3in x .5in
- 750-630 LEEP / LLETZ Sm 7.5cm x 2cm/3in x .75in
- 750-631 LEEP / LLETZ Sm W/Smoke Tube 7.5cm x 2cm/3in x .75in
- 750-632 LEEP / LLETZ Sm W/DSE Tube 7.5cm x 2cm/3in x .75in
- 750-633 LEEP / LLETZ Med 10cm x 4cm /4in X 1.5in
- 750-634 LEEP / LLETZ Med W/Smoke Tube 10cm x 4cm /4in X 1.5in
- 750-635 LEEP / LLETZ Med W/DSE Tube 10cm x 4cm /4in X 1.5in
- 750-636 LEEP / LLETZ Lrg 12cm x 4cm/4.75in x 1.5in
- 750-637 LEEP / LLETZ Lrg W/Smoke Tube 12cm x 4cm/4.75in x 1.5in
- 750-638 LEEP / LLETZ Lrg W/DSE Tube 12cm x4cm/4.75in x 1.5in

## PEDERSON SPECULUM

- 750-639 LEEP / LLETZ Sm 9cmx1.5cm/3.5in x .5in
- 750-640 LEEP / LLETZ Sm W/Smoke Tube 9cmx1.5cm/3.5in x .5in
- 750-641 LEEP / LLETZ Sm W/DSE Tube 9cmx1.5cm/3.5in x .5in
- 750-642 LEEP / LLETZ Med 10cm x 2.5cm/4in x 1in
- 750-643 LEEP / LLETZ Med W/Smoke Tube 10cm x 2.5cm/4in x 1in
- 750-644 LEEP / LLETZ Med W/DSE Tube 10cm x 2.5cm/4in x 1in
- 750-645 LEEP / LLETZ Lrg 11.5cm x 2.5cm/ 4.5in x 1in
- 750-646 LEEP / LLETZ Lrg W/Smoke Tube 11.5cm x 2.5cm/ 4.5in x 1in
- 750-647 LEEP / LLETZ Lrg W/DSE Tube 11.5cm x 2.5cm/ 4.5in x 1in

## VU-MAR SPECULUM (4CM YOKE)

- 750-651 LEEP / LLETZ Sm 9cm x 2.5cm/3.5in x 1.0in
- 750-652 LEEP / LLETZ Med W/Smoke Lrg Tube 10cmx2.5cm/4in x 1in
- 750-653 LEEP / LLETZ Lrg W/Smoke Tube 12cm x 4cm/4.75in x 1.5in

## VU-MAR XL SPECULUM GRAVES (SQUARE YOKE)

- 750-654 LEEP / LLETZ Vu-Mar XL Speculum Graves Med W/Smoke Tube 10cm x 4cm /4.0in x 1.5in
- 750-655 LEEP / LLETZ Vu-Mar XL Speculum Graves Lrg W/Smoke Tube 12cm x 4cm/4.75in x 1.5in





# SPECULA

A wide array of Elektrobue Leep Speculums in many different shapes and sizes allow for easy access and viewing of the surgical site, available w/your choice of smoke evacuation options.

## **VU-MAR XL SPECULUM PEDERSON**

750-656 LEEP / LLETZ Med W/Smoke Tube 10cm x 2.5cm/4.0in x 1.0in

750-657 LEEP / LLETZ Lrg W/Smoke Tube 12cm x 2.5cm/4.75in x 1.0in

## **WEISMAN GRAVES SPECULUM**

750-658 LEEP / LLETZ Left Opening Med W/Smoke Tube 10cm x 4cm/4in x 1.5in

750-659 LEEP / LLETZ Left Opening Lrg W/Smoke Tube 12cm x 4cm/4.75in x 1.5in

750-660 LEEP / LLETZ Right Opening Med W/Smoke Tube 10cm x 4cm/4in x 1.5in

750-661 LEEP / LLETZ Right Opening Lrg W/Smoke Tube 12cm x 4cm/4.75in x 1.5in

## **WEISMAN PEDERSON SPECULUM**

750-662 LEEP / LLETZ Right Opening Med 10cmx2.5/4inx1.0in

750-663 LEEP / LLETZ Right Opening Med W/Smoke Tube 10cmx2.5/4inx1.0in

750-664 LEEP / LLETZ Right Opening Med W/DSE Tube 10cmx2.5/4inx1.0in

750-665 LEEP / LLETZ Righth Opening Lrg 12cmx2.5cm/4.75inx 1in

750-666 LEEP / LLETZ Righth Opening Lrg W/Smoke Tube 12cmx2.5cm/4.75inx 1in

750-667 LEEP / LLETZ Righth Opening Lrg W/DSE Tube 12cmx2.5cm/4.75inx 1in

## **SEGOND SPECULUM**

750-668 LEEP / LLETZ 12cm/ 4.75in

## **STEINER SPECULUM**

750-669 LEEP / LLETZ 1.1 Kg/2.5 lb 14cm x 3cm/5.5in x 1.25in

## **TAUBER SPECULUM**

750-670 LEEP / LLETZ 9cm/3 1/2in

# KIT

This customizable set of instruments used for colposcopy comes with all the instrumentation needed for any practice.



## **KIT**

KIT-001 Nurse Practitioner Mobile LEETZ Kit

### **CHOOSE ONE OF EACH INSTRUMENT**

Biopsy Punch Tischler/Kevorkian

Speculum Graves/Pederson

Endocervical speculum W/Ratchet/W/Set Screw

Hooks Iris/Emmett

Endocervical currette W/Basket/Without Basket



# RETRACTORS

Elektroblue Leep Retractors are the best choice to help provide a clearer visual of the cervix.

## **BREISKY RETRACTOR**

750-701 100mm x 18mm/ 4in x .70in

750-702 100mm x 25mm/4in x 1in

750-703 100mm x 30mm/4in x 1.25in

750-704 100mm x 35mm/4in x 1.3in

750-705 100mm x 38mm/4in x 1.5in

750-706 120mm x 30mm/4.72in x 1.25in

750-707 130mm x 35mm/ 5in x 1.3in

750-708 140mm x 35mm/ 5.5in x 1.3in

750-709 150mm x 35mm/6in x 1.3in

750-710 160mm x 35mm/6.3in x 1.3in

## **BREISKY-NAVRATIL RETRACTOR**

750-711 170mm x 25mm / 6.75in x 1in

## **DEAVER RETRACTOR**

750-712 Size 25mm x 23cm/1in x 9in

750-713 Size 25mm x 30cm/1in x 12in

## **DESMARRES LID RETRACTOR**

750-714 10mm 14cm/5.5in

750-715 12mm 14cm/5.5in

750-716 14mm 14cm/5.5in

750-717 16mm 14cm/5.5in

750-718 Double Ended 8/10mm 15cm /6in

## **FARABEUF RETRACTOR**

750-719 Set of Two 15cm/6in

## **FOMON RETRACTOR**

750-720 Double Ball End 16cm/6.25in

## **GELPI RETRACTOR**

750-721 11.5cm/4.5in

## **JACKSON RETRACTOR**

750-722 SM 3.8cm x 7.6cm /1.5in x 3in

750-723 MED Single Ended 3.8cm x 8.9cm/1.5in x 3.5in

750-724 LG Single Ended 3.8cm x 10.2cm/1.5in x 4in

## **JAEGER LID RETRACTOR**

750-725 11cm/4.25in

## **LATERAL WALL RETRACTOR**

750-726 Long & Narrow 5.75cm/2.25in Blades

750-727 Lateral Wall Retractor 7.5cm/2.5in Blades

750-728 Lateral Wall Retractor 8.125in

## **RAGNELL DAVIS RETRACTOR**

750-729 15cm/6in

## **REDUNDANT WALL RETRACTOR**

750-730 9cm x 2cm /3.5in x .75in Blades 18cm/7in

750-731 12.5cm X 2cm /5in X .75in Blades 18cm/7in

## **RIBBON RETRACTOR**

750-732 Malleable 6mm: 18cm/7in

750-733 Malleable 10mm: 18cm/7in

## **RIGBY RETRACTOR**

750-734 .75in X 3.5in Blades 7.75in

## **SENN MUELLER RETRACTOR**

750-735 Sharp 16cm/6.25in

750-736 Blunt 16cm/6.25in

## **SIMS RETRACTOR**

750-737 MED Single Ended

750-738 MED DSE Smoke Tube Single Ended

750-739 MED Double Ended

750-740 MED DSE Smoke Tube Double Ended

## **U.S. ARMY RETRACTOR**

750-741 Parker Langenbeck Set of Two 21.5cm/8.5in

