



Precision. Performance. Peace of Mind.

SurgyCut

Fully Disposable Laparoscopic Graspers



Sold as all-in-one pre-assembled device, including handle, shaft, and tip. Ratcheted and non-ratcheted handle options available!

Monopolar Energy | Sterile-Packed

5mm Shaft | 35 & 45 cm lengths

- **Single-Use Simplicity** - No reprocessing, no cross-contamination, no downtime. Every case starts with an all-in-one (handle, shaft, & tip) clean, sharp instrument.
- **Surgeon-Centered Ergonomics** - Lightweight, balanced handle with smooth actuation and tactile feedback—designed for long cases and fine control.
- **Secure Tissue Handling** - High-precision jaws provide consistent grip strength and anti-slip patterns.

All graspers sold in boxes of 10 ea

Add "-R" to end of item number for Ratchet Handle version.

Item Number	Jaw Pattern	Image
790-354-535 (35 cm)	Allis	
790-354-545 (45 cm)		
790-352-535 (35 cm)	Babcock	
790-352-545 (45 cm)		
790-355-535 (35 cm)	Bowel	
790-355-545 (45 cm)		
790-353-535 (35 cm)	Claw	
790-353-545 (45 cm)		
790-351-535 (35 cm)	Clinch	
790-351-545 (45 cm)		
790-350-535 (35 cm)	Fenestrated	
790-350-545 (45 cm)		
790-300-535 (35 cm)	Maryland	
790-300-545 (45 cm)		
790-302-535 (35 cm)	Mixer	
790-302-545 (45 cm)		